

MSDS Red Cell Reagents in CellStab

Material Safety Data Sheet

Red Cell Reagents in CellStab

1. Identification of the substance/preparation and company

Human red cell reagents for *in-vitro* diagnostic use only:

ID Panel	PR142 PR143
ID Panel enzyme	PR152 PR153
A ₁ rr cell	PR014
Brr cell	PR035
2 cell screen	PR102
3 cell screen	PR122
Titration Cells	PR407
2 cell rr screen	PR106

NBS Reagents

14 Estuary Banks, Speke, Liverpool L24 8RB

Telephone: 0151 552 7157

Fax: 0151 552 7156

2. Composition/information on ingredients

Donated human red cells as a 0.8% cell suspension in CellStab which is Glycine buffered saline containing sugars, trimethoprim and sulfamethaxazol as preservatives as supplied by DiaMed.

3. Hazards identification

CellStab solution contains no known hazardous material.

All blood grouping reagents should be treated as potentially infectious. The donations used in these products are not sterilised; capable of transmitting any biological agent that has not been detected by routine screening at the time of manufacture. No known test methods can offer assurances that products derived from human blood will not transmit infectious disease. Appropriate care should be taken in the use and disposal of this product.

4. First aid measures

Eye contamination: Immediate and prolonged irrigation with copious amounts of water. Seek medical advice

Skin contamination: Wash thoroughly with copious amounts of water. If skin is broken seek medical advice

Ingestion: Wash mouth with copious amounts of water and seek medical advice

Inhalation: N/A

5. Fire fighting measures

Non combustible. If involved in fire use extinguishing media appropriate to the surrounding conditions.

MSDS Red Cell Reagents in CellStab

- 6. Accidental release measures**
Wear laboratory coats and disposable nitril gloves. Control spillage in accordance with HSE regulations. Spray spillage with appropriate virucidal detergent and absorb, wash spill site after material pick up is complete.
- 7. Handling and storage**
When handling product wear laboratory coats and disposable nitril gloves. Store the reagent at +2°C to +8°C.
Protect from contamination.
- 8. Exposure controls/personal protection**
When handling product wear laboratory coats and disposable nitril gloves.
- 9. Physical and chemical properties**
A suspension of red cells in a clear liquid.
- 10. Stability and reactivity**
Do not freeze
Protect from extremes of temperature
Do not use if reagent exhibits obvious discolouration or haemolysis.
- 11. Toxicological information**
Data not available
- 12. Ecological Information**
Data not available
- 13. Disposal Considerations**
Dispose of in accordance with Hazardous Waste Regulations 2005 and any additional local regulations
Also see 3
- 14. Transport Information**
Suitable for transportation by road and rail.
Packaging complies with regulations in CHIP (The Chemicals (Hazard Information and Packaging for Supply) Regulations 2002).
- 15. Regulatory information**
For use *in-vitro* diagnostic tests by appropriately trained laboratory personnel.
COSHH
Health and Safety at Work Act 1974
The Chemicals (Hazard Information and Packaging for Supply) Regulations 2002
Environmental Protection Act 1990
Hazardous Waste Regulations 2005
EEC/91/155
- 16. Other information**
The information in this safety data sheet does not replace the users own assessment of work place risk as required by other health and safety legislation.