

ID Panel Profile Products PR141, 142, 143, 144, 146, 151, 152, 153 & 154



IVD		NBS REAGENTS		PANEL 1	
Product	Lot No.	Product	Lot No.	Product	Lot No.
ID Panel in Alsevers	R144 3172	ID Panel Papainised in Alsevers	R154 3172	ID Panel in LISP	R146 3172
ID Panel in CellStab	R143 3172	ID Panel Papainised in CellStab	R153 3172	Expiry Date : 2009.07.30	

Patient's Name	Ref. No.	Sample No.	Conclusion
D.O.B.			Tested by: _____ Date: _____

Unless otherwise indicated, all cells are positive for Kp^b and Lu^b and negative for Wr^a and Co^b.

Instructions for use can be found at http://www.blood.co.uk/hospitals/diagnostic_services/reagents/index.asp#Pro

	Rh	C	D	E	c	e	C ^w	M	N	S	s	P1	Lu ^a	K	k	Kp ^a	Le ^a	Le ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Other															
1	R ₁ ^w R ₁	+	+	0	0	+	+	0	+	0	+	1	+	0	+	0	0	+	+	0	+	0																
2	R ₁ R ₁	+	+	0	0	+	0	+	0	+	0	3	0	+	+	0	+	0	0	+	0	+																
3	R ₂ R ₂	0	+	+	+	0	0	0	+	0	+	0	0	0	+	0	+	0	+	0	0	+																
4	r'r	+	0	0	+	+	0	+	+	+	0	2	0	0	+	0	0	+	+	0	+	0																
5	r''r	0	0	+	+	+	0	+	0	+	0	1	0	0	+	0	0	+	0	+	+	0																
6	rr	0	0	0	+	+	0	+	+	0	+	4	0	+	0	0	0	0	0	+	+	0																
7	rr	0	0	0	+	+	0	0	+	0	+	1	0	+	+	0	+	0	+	0	0	+																
8	rr	0	0	0	+	+	0	+	0	+	0	2	0	0	+	+	+	0	0	+	0	+																
9	rr	0	0	0	+	+	0	0	+	+	0	3	0	0	+	0	0	+	0	+	+	+																
10	rr	0	0	0	+	+	0	+	+	0	+	0	+	0	+	0	0	+	+	0	+	0																
Reagent																																						
Lot No.																																						