

3 Cell Screen Profile Product PR121 & PR122



IVD

NBS REAGENTS LIVERPOOL

Product	Lot No	Product	Lot No.	Expiry Date
Alsevers	R121 3281	CellStab	R122 3281	2009.07.16

Unless otherwise indicated, all cells are positive for Kp^b and Lu^b and negative for Wr^a Lu^a and Co^b

Instructions for use can be found at http://hospital.blood.co.uk/hospitals/diagnostic_services/reagents/index.asp#Pro

	Rh	C	D	E	c	e	C ^w	M	N	S	s	P1	K	k	Kp ^a	Le ^a	Le ^b	Fy ^a	Fy ^b	JK ^a	JK ^b	Other
1	R₁^wR₁	+	+	0	0	+	+	+	+	0	+	+	+	+	0	0	+	+	0	+	0	Cob+
2	R₂R₂	0	+	+	+	0	0	+	0	0	+	0	0	+	0	0	+	0	+	0	+	
3	rr	0	0	0	+	+	0	+	0	+	0	+	0	+	+	+	0	+	0	+	+	