

## Papainised 3 Cell Screen Profile



IVD



### NHSBT REAGENTS

Product	Lot No	Expiry Date
<b>Papainised in Alsevers</b>	<b>R124 3530</b>	<b>2018.04.19</b>

Instructions for use can be found at [www.blood.co.uk/reagents](http://www.blood.co.uk/reagents)

	Rh	C	D	E	c	e	C <sup>w</sup>	P1	K	k	Le <sup>a</sup>	Le <sup>b</sup>	Jk <sup>a</sup>	Jk <sup>b</sup>	Other
<b>1</b>	<b>R<sub>1</sub><sup>w</sup>R<sub>1</sub></b>	<b>+</b>	<b>+</b>	<b>0</b>	<b>0</b>	<b>+</b>	<b>+</b>	<b>w</b>	<b>0</b>	<b>+</b>	<b>0</b>	<b>+</b>	<b>+</b>	<b>+</b>	
<b>2</b>	<b>R<sub>2</sub>R<sub>2</sub></b>	<b>0</b>	<b>+</b>	<b>+</b>	<b>+</b>	<b>0</b>	<b>0</b>	<b>+</b>	<b>+</b>	<b>+</b>	<b>0</b>	<b>+</b>	<b>+</b>	<b>0</b>	
<b>3</b>	<b>rr</b>	<b>0</b>	<b>0</b>	<b>0</b>	<b>+</b>	<b>+</b>	<b>0</b>	<b>+</b>	<b>0</b>	<b>+</b>	<b>+</b>	<b>0</b>	<b>0</b>	<b>+</b>	