





Component Description	Shelf life	NHSBT Pulse Code	Barcode	Component Status
Red Cells, Washed LD Manual Preparation Suspended in SAGM	14 days	0076	 a0465313b	<b>NEW</b>  Please add to LIMS system
Red Cells, Washed, LD, Irradiated Manual Preparation Suspended in SAGM	Max 48 hours	G076	 a0465323b	<b>NEW</b>  Please add to LIMS system
Red Cells, Washed, LD Automated Preparation	14 days	0074	 a0100103b	Discontinued.
Red Cells, Washed, LD, Irradiated Automated Preparation	Max 48 hours	G074	 a0200103b	Discontinued.